



SemenHos

FOR PROFESSIONAL USE ONLY



Application

This SemenHos test is used to test the vitality of sperm cells. The hypo-osmotic swelling is based on the semi-permeability of the intact cell membrane and their ability of active water transport, in order not to burst. In sperm with intact membranes the flagellum swells up within 5 min. This change remains stable for up to 30 min.

Principle

In this hypo-osmotic swelling test swelling of cells only occurs in vital cells with an intact membrane by using hypotonic solution.

Storage and stability

2-8 °C. Sterile sampling. Contains no antibiotic.

24 months from date of manufacture.
After opening use within 7 days.

Content

SemenHos solution, 20 x 900 µl

Necessary utensils

- Coverslips (18 x 18 mm)
- Gloves
- Contrasting phase microscope
- Native ejaculate or washed sperm (105-110 µl)
- Slides
- Paper towels
- Pipettes and tips (10-100 µl)
- Water bath or heating cabinet (37°C)

Preparation of SemenHos solution

Preheat the SemenHos solution to 37 °C

Procedure

1. Occasionally native ejaculate without additions includes hypo-osmotic sperm forms. Transfer 5-10 µl liquefied semen without air bubbles to a slide and cover it with cover slips. Microscope at 400x magnification. This is the zero-value.
2. Examine the percentage of sperm with swollen flagellum by observing 100 sperm, calculated in duplicate. Note this value note* (a_0 %).
3. Add 100 µl ejaculate without air bubbles to 900 µl preheated SemenHos solution and mix
4. Incubate this mixture for 10 min at 37 °C.
5. Transfer 10 µl of the mixture to a slide and place and cover it with a cover slip.
6. Microscope at 200x or 400x magnification.
7. Repeat twice step 2 to 6.




Fig.1: Schematic representation of typical morphological changes of human spermatozoa after exposure with hypo-osmotic solution (extract WHO 2010).

No change (a). Different tail changes (b-g).
The swelling in the tail region is indicated in gray.

Gi183/V1

Planer Limited

110 Windmill Road, Sunbury-On-Thames
Middlesex TW16 7HD, United Kingdom
A Hamilton Thorne Company

Tel: +44 (0)1932 755 000
enquiries@planer.com
www.planer.com

Manufactured by

ANTIGENES GmbH

Labordiagnostika
Römerstraße 20-22
D-58332 Schwelm, Germany

Tel.: +49 2336 9154958
info@antigenes.de
www.antigenes.de



Evaluation

Percentage of vitality of spermatozoa

Calculate the difference between the percentage of sperm with swollen flagellum before and after incubation with the SemenHos solution.

Example:

Before incubation:

1. Count 2/100

2. Count 3/100

Mean value: $2.5 / 100 = 2.5\%$

After incubation:

1. Count 88/100

2. Count 97/100

Mean value: $92.5 / 100 = 92.5\%$

Result: $92.5\% - 2.5\% = 90\%$

90 % of the sperm are vital

The SemenHos test is regarded as normal, if after incubation more than 60% of the sperm show a swollen flagellum.

The sample is not normal when the result is less than 50% (WHO 2010).

Safety information / Precautions

(Please read also safety data sheets)

- All semen samples should be considered potentially infectious. Handle all samples as though they are HIV or hepatitis infected material.
- When working with samples and reagents always wear protective clothing (gloves, gowns, eye /face protection).
- All ingredients of reagents are classified as non-toxic

References

1. Drevious L, Eriksson H, (1966) Osmotic swelling of mammalian spermatozoa, *Experimental Cell Research*, 42: 136-56
2. Jeyendran RS, et al, (1984) Development of an assay to assess the functional integrity of the human sperm membrane and its relationship to the other sperm characteristics, *Journal of Reproduction and Fertility*, 70: 219-28
3. WHO Press, (2010) Laboratory manual for the examination and processing of human semen, 5th edition
4. Zaneveld LJD, (1984) *Journal of Reproduction and Fertility*, 70:219-228. © Society for Reproduction and Fertility

Planer Limited

110 Windmill Road, Sunbury-On-Thames
Middlesex TW16 7HD, United Kingdom
A Hamilton Thorne Company

Tel: +44 (0)1932 755 000
enquiries@planer.com
www.planer.com

Manufactured by ANTIGENES GmbH

Labordiagnostika
Römerstraße 20-22
D-58332 Schwelm, Germany

Tel.: +49 2336 9154958
info@antigenes.de
www.antigenes.de