



Vitrification




Preparation

Ensure all media are warmed up to 37 °C and mixed well before use.

We would strongly advise you to read through all the steps of the vitrification/thawing procedure before starting the procedure.

Thawing

- In a 4-well dish fill the first well with 300 µl of VitriStore ThawMedium 1, the second with VitriStore ThawMedium 2, the third with VitriStore ThawMedium 3 and the 4th with VitriStore ThawMedium 4.



- Remove the vitrification straw from the outer sheath as indicated in the instructions for use of the vitrification device.

- Immediately plunge the vitrification straw into pre-heated VitriStore Thaw Medium 1 (37 °C) and leave in Thawing 1 medium for 3 minutes.



- Transfer into VitriStore Thaw Medium 2 (37 °C) and leave in this medium for 2 minutes.

- Transfer into VitriStore Thaw Medium 3 (37 °C) and leave in this medium for 2 minutes.

- Finally transfer into VitriStore Thaw Medium 4 (37 °C) and wash for at least 1 minute.

- Transfer into blastocyst culture medium for continued cell culture (e.g. GM501 Cult media).

35 mm Petri dish	Well 1	Well 2	Well 3
VST 1	VST 2	VST 3	VST 4
3 min.	2 min.	2 min.	1 min.

Gi136/V1

Planer Limited

110 Windmill Road, Sunbury-On-Thames
Middlesex TW16 7HD, United Kingdom
A Hamilton Thorne Company

Tel: +44 (0)1932 755 000
enquiries@planer.com
www.planer.com

Manufactured by Gynemed GmbH & Co. KG

Lübecker Straße 9 23738 Lensahn
Germany

Tel.: +49 (0) 4363 90 32 90
info@gynemed.de
www.gynemed.de