
# Vitrification Cryoprotocol for Oocyte and Embryo

Code:VT601, VT602

**EN: Intended use**  
This product is to be used for vitrification of oocytes (MI) and embryos.

**Vitrification Media**  
No. 1 (Color code: white) Basic Solution (BS): 1x1.5 ml vial only for oocyte (MI) vitrification.  
No. 1 (Color code: blue) Equilibration Solution (ES): 1x1.5 ml vial only for oocyte (MI) vitrification.  
No. 2 (Color code: green) Vitrification Solution (VS): 2x1.5 ml vials only for oocyte (MI) vitrification.  
No. 1 (Color code: white) Basic Solution (BS): 1x1.5 ml vial only for embryo (MI) vitrification.  
No. 1 (Color code: blue) Equilibration Solution (ES): 1x1.5 ml vial only for embryo (MI) vitrification.  
No. 2 (Color code: green) Vitrification Solution (VS): 2x1.5 ml vials only for embryo (MI) vitrification.

**Instructions for use (IFU)**  
Preparation  
Fill 90% of the Cooling Rack with fresh liquid nitrogen.  
Use the rack with the thickness of zona pellucida and record it.  
Bring BS, ES and VS to room temperature (23-27°C as recommended).  
Use certified pipette as a handling tool, with a suitable internal diameter for oocytes or embryos.  
The recommended internal diameters are as follows: 120 µm for oocytes (MI), 120 µm for pronuclear stage embryos, 150 to 180 µm for cleavage stage embryos, and 180 to 250 µm for blastocyst stage embryos.

**Oocyte**  
1. Equilibration procedures for Oocyte (MI) and Embryo are different.  
2. Equilibration of oocyte (MI)  
1. Drop 20 µl of BS into first well and each 300µl of VS into second and third well on the Repro Plate with pipette.  
2. Step 0: Transfer the Embryo to the TOP center of the BOTTOM of BS.  
3. Step 1: Immediately add ES 20 µl gently on the TOP of BS of the first well and leave it for 3 minutes.  
4. Add another ES 20 µl gently on the TOP of the first well and leave it for 3 minutes.  
5. Add another ES 240 µl gently on the TOP of the first well and leave it for 9 minutes.

**Embryo**  
1. Drop each 300 µl of ES into first well, VS into second and third well on the Repro plate with pipette.  
2. Step 0: Transfer the Embryo to the TOP center of BS of the culture dish. It will spontaneously begin to shrink and then gradually return to its original size by absorbing the ES solution (within 15 minutes).

**Quality control specification**  
The following tests were performed for each lot of this product:  
Sterility by Sterility Test (EP)  
Endotoxin by LAL method.  
Mouse Embryo Assay (One Cell) pH (EP)  
Osmolality (EP)

**Storage instructions and stability**  
Store the vials at 2 to 8 °C.  
This product is stable until the expiry date labeled on the vial.

**Composition**  
-HEPES with Basic Culture Medium  
-Ethylene Glycol  
-Dimethylsulfoxide  
-Trehalose  
-Hydroxypropyl cellulose  
-Gentamicin

**Contraindications**  
Do not re-sterilize.  
Do not use. Re-use may cause a contamination.  
Do not use solution that shows cloudiness or becomes decomposed.  
Do not use the product if you notice any unusual regarding the solution on the label (number, color, name, volume).  
Morphologically abnormal oocytes, embryos, or significantly poor grade oocytes are unsuitable for cryopreservation.  
Observe all federal, state and local environmental regulations when discarding the product.

**Warning**  
Read the instructions for use prior to use.  
Use the KITAZATO vitrification media to thaw the oocytes (MI) or embryos.  
Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.  
NOTE: The long-term safety of a cryoprotection technique and maximum storage in liquid nitrogen has not been established and unknown.

Copyright © Kitazato Corporation All Rights Reserved.

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