



# SemenStain

FOR PROFESSIONAL USE ONLY



## Application

SemenStain is a quick-staining-method to assess the morphology of sperm (spermioqram). This method is composed of a staining kit which allows differential staining of the sperm parts due to their different basophilic, eosinophilic and neutrophilic properties.

## Principle

The sperm are fixed. Here, the succedan-staining is used, which means, three dyes are used one after the other and it results to differentiate staining of different tissues with individual dyes.

## Storage and stability

15-25 °C

36 months from date of manufacture

## Content

- Reagent 1, 1x 50 or 250 ml
- Reagent 2, 1x 50 or 250 ml
- Reagent 3, 1x 50 or 250 ml
- Reagent 4, 1x 50 or 250 ml

## Necessary utensils (not included)

- Native ejaculate or washed sperm (5-10 µl)
- Staining cuvettes (8x) or tubes (50 ml, 8x)
- Gloves
- Tweezers
- Paper towels
- Slides
- Slides rack (if more than five slides are to dye)
- Immersion oil
- Microscope

## Procedure (see also diagram)

1. Apply 5-10 µl sperm per slide, smear the sperm with a coverslip and let dry. We recommend preparing 2 slides per patient.
2. Fill the staining cuvettes with reagent 1, reagent 2, reagent 3 and reagent 4, respectively. Fill 4 other empty cuvettes with water. Place the cuvettes side by side. Label them from 1 to 8.
3. Immerse slides 3 minutes by repeated immersion in cuvette 1 (reagent 1) to fix the preparation. Wash slides 3 minutes in cuvette 2 (water). Then place the slides vertically on paper towels to remove excess water.
4. Colour the slides 1 minute by repeated immersion in cuvette 3 (reagent 2). Wash the slides in cuvette 4 (water). Change the water several times until water stays clear; remove with paper towel the excess water of the slides.
5. Colour the slides 1 minute by repeated immersion in cuvette 5 (reagent 3), wash the slides according to step 4 in cuvette 6 (water) and remove from the slides with paper towel the excess water.
6. Colour the slides 1 minute by repeated immersion in cuvette 7 (reagent 4), wash the slides according to step 4 in cuvette 8 (water) and remove from the slides with paper towel the excess water.
7. Dry the slides in the open air.

## Evaluation

Evaluate the sperm with immersion oil at 1000x magnification on the side of the slide with lower sperm density. Here the sperm are better to assess individually.

Gi186/V1

### Planer Limited

110 Windmill Road, Sunbury-On-Thames  
Middlesex TW16 7HD, United Kingdom  
A Hamilton Thorne Company

Tel: +44 (0)1932 755 000  
enquiries@planer.com  
www.planer.com

### Manufactured by ANTIGENES GmbH

Labordiagnostika  
Römerstraße 20-22  
D-58332 Schwelm, Germany

Tel.: +49 2336 9154958  
info@antigenes.de  
www.antigenes.de



The criteria for classification of sperm by their morphology can be found in the WHO laboratory manual (2010).

Sperm cell parts	Staining/colour
Head - Nucleus	red
Head - Acrosome	dark green
Middle part	pale green
Tail	green

After the evaluation the immersion oil can be gently removed from the slide with a paper towel. Then the slide can be immersed in reagent 1 for 5 min, dried and stored. It is also possible to produce preparations with a coverslip and glue for long term storage.


**References**

- Andersen AG et al. (2000) High frequency of sub-optimal semen quality in an unselected population of young men. *Human Reproduction*, 15:366-372
- Auger J, Eustache F (2000) Standardisation de la classification morphologique des spermatozoïdes humains selon la méthode
- Behre HM et al. (2000) Diagnosis of male infertility and hypogonadism. In: Nieschlag E, Behre HM, eds. *Andrology, male reproductive health and dysfunction*. Springer: 92ff.
- Cooper TG et al. (2010) World Health Organization reference values for human semen characteristics. *Human Reproduction Update*, 16:231-245
- Cross NL (1995) Methods for evaluating the acrosomal status of human sperm. In: Fenichel P, Parinaud J, eds. *Human sperm acrosome reaction*. Paris, John Libbey Eurotext (Colloques INSERM): 277-285
- Kruger TF et al. (1987) A quick, reliable staining technique for human sperm morphology. *Archives of Andrology*, 18:275-277
- WHO Press, (2010) *Laboratory manual for the examination and processing of human semen*. 5th edition

**Safety information / precautions**

- All semen samples should be considered potentially infectious.
- Handle all samples as though they are HIV or hepatitis infected material.
- When working with samples and reagents always wear protective clothing (gloves, gowns, eye / face protection).
- Reagent 1 is containing methanol: toxic by inhalation, skin contact or ingestion. May cause organ damage. There is a risk of irreversible damage.
- All other ingredients are not classified as toxic

**Diagram (see also procedures outlined above)**



Gi186/V1

**Planer Limited**  
 110 Windmill Road, Sunbury-On-Thames  
 Middlesex TW16 7HD, United Kingdom  
 A Hamilton Thorne Company

**Tel:** +44 (0)1932 755 000  
 enquiries@planer.com  
 www.planer.com

**Manufactured by**  
**ANTIGENES GmbH**  
 Labordiagnostika  
 Römerstraße 20-22  
 D-58332 Schwelm, Germany

Tel.: +49 2336 9154958  
 info@antigenes.de  
 www.antigenes.de