

**KITAZATO**

IVF Needles

Oocyte Retrieval Needles

Needles designed to maximise control and patient comfort.

The Kitazato range of oocyte retrieval needles, with their steel triple cut blades, are designed to provide a fast puncture, minimising damage to the ovarian tissue.

Good control of the puncture, allows a reduced operation time helping oocytes to be retrieved in the best possible condition.

Kitazato OPU needles offer the widest range of diameters on the market – from 16 G to 21 G. This provides ideal options for working with patients with low response in natural cycles and allows follicular aspiration to be undertaken, without an anaesthetic, using the smallest diameters.

- Improved patient comfort
- Maximum control
- Reduced procedure time
- Large range of sizes
- Single or dual lumen



Oocyte Retrieval Needles

Kitazato OPU needles offer the widest range of diameters on the market from 16 G to 21 G

| Order Number | Description | Guage | Length | Aspiration Line |
|----------------------------------|-----------------------------|-------|---------|---------------------------|
| Pump Aspiration Single Lumen | | | | |
| FDM326350 | | 16 G | 35 cm | 80 cm* |
| FDM327350 | | 17 G | 35 cm | 80 cm* |
| FDM328350 | | 18 G | 35 cm | 80 cm* |
| FDM329350 | | 19 G | 35 cm | 80 cm* |
| FDM320300 | | 20 G | 30 cm | 80 cm* |
| FDM320350 | | 20 G | 35 cm | 80 cm* |
| FDM321300 | | 21 G | 30 cm | 80 cm* |
| Manual Aspiration Single Lumen | | | | |
| FDM316001 | Manual Aspiration Luer Lock | 16 G | 32.5 cm | |
| FDM317001 | Manual Aspiration Luer Lock | 17 G | 32.5 cm | |
| FDM318001 | Manual Aspiration Luer Lock | 18 G | 32.5 cm | |
| FDM319001 | Manual Aspiration Luer Lock | 19 G | 32.5 cm | |
| FDM310001 | Manual Aspiration Luer Lock | 20 G | 32.5 cm | |
| FDM367350 | Manual Aspiration Syringe | 17 G | 35 cm | 60 cm* |
| Kizato OPU Needle - Double Lumen | | | | |
| FDM376350 | Double Lumen | 16 G | 35 cm | 100 cm + 80 cm (flushing) |
| FDM378350 | Double Lumen | 17 G | 35 cm | 100 cm + 80 cm (flushing) |
| FDM377350 | Double Lumen | 18 G | 35 cm | 100 cm + 80 cm (flushing) |
| Other | | | | |
| FDM357353 | 3 way stopcock | 17 G | 35 cm | 80 cm + 40 cm* |
| FDM347275 | IVM needle | 17 G | 27.5 cm | 80 cm |

* Different lengths of aspiration line available. ** Kitazato OPU needles are sold in boxes of 10 units.

Quality Control

- MEA (Mouse Embryo Assay): One cell assay \geq 80% after 96 hours
- Endotoxin: < 20EU/device
- SAL 10^{-6}
- Cytotoxicity Test
- Intracutaneous Reactivity Test
- Sensitisation Text

Composition

- Stainless Steel SUS304

Ki143/V1

Planer Limited

110 Windmill Road, Sunbury-On-Thames
Middlesex TW16 7HD, United Kingdom
A Hamilton Thorne Company

Tel: +44 (0)1932 755 000
enquiries@planer.com
www.planer.com

Manufactured by Kitazato Corporation

81 Nakajima, Fuji
Shizuoka 416-0907

Tel: +81-3-3434-2731
trading@kitazato.co.jp
www.kitazato.co.jp