



Cryogenic RFID Witnessing System

SINGLE READER MK2

THE SINGLE READER MK2 IS THE NEXT ITERATION OF THE SINGLE READER SERIES.

It now uses a USB C port for the power supply and data transfer. It also incorporates the latest RFID reader module.

SPECIFICATIONS

Communication	USB C
Power Requirement	Power via USB
Dimensions L x W x H (mm)	117 x 100 x 32
Weight	92g
LED	Green LED indicated if unit is powered on
RFID	ISO15693 High Frequency 13.56 MHz I-CODE1 and ISO/IEC 15693 protocols
Data Rates	115.2 kbit/s
Protection Rating	IP 54
Operating Temperature	-40 to +80 °C
Storage Temperature	0 to +35 °C
Case Material	Nylon 12
Software	Drivers compatible with Windows 10 and 11
Warranty	1 year

WANT TO FIND OUT MORE ABOUT CRYOGENIC RFID? CONTACT US.

Specifications may change without notice. Third party trademarks acknowledged.

Si213/V1

© 2023 Planer Ltd

PLANER
PRESERVE / PROTECT / NURTURE

A Hamilton Thorne Company

www.planer.com

Manufactured by:
Cryogatt Systems Ltd
7, Abbots Business Park
Primrose Hill, Kings Langley
WD4 8FR, United Kingdom

Planer Limited
110 Windmill Road
Sunbury-On-Thames
TW16 7HD, United Kingdom
+44 (0)1932 755 000

Contact us: planer.com/contact