



## Micromanipulators

# MN-4VM Manual Coarse Manipulator

### Manual coarse manipulator suitable for Setagaya

This manipulator is for the quick positioning of the pipette into the microscope field of view.

Without the need for electricity, this setup is more economical than its motorised counterpart.

- Quick positioning
- Responsive operation
- Performance proven design
- Different combinations for a variety of installations



# MN-4VM Manual Coarse Manipulator for Setagaya

This manipulator is for the quick positioning of the pipette into the microscope field of view



Order code	
FNHMN-4	MN-4VM Setagaya manual coarse manipulator

Specification	
Accessories included	Hex Wrench
Movement range	Coarse X30 mm, Y30 mm, Z30 mm
Dimensions and Weight	24 (W) × 55 (D) × 54 (H) mm, 58 g

Ni192/V1

## Planer Limited

110 Windmill Road, Sunbury-On-Thames  
Middlesex TW16 7HD, United Kingdom  
A Hamilton Thorne Company

**Tel:** +44 (0)1932 755 000  
enquiries@planer.com  
www.planer.com

## Manufactured by

Narishige Scientific Instrument Lab.  
27-9, Minamikarasuyama 4-chome,  
Setagaya-ku, Tokyo 157-0062  
Japan

**Tel:** +81-(0)3-3308-8233  
sales@narishige.co.jp  
www.narishige-group.com