



Cryopreservation

GM501 EmbryoStore

GM501 EmbryoStore is a set of ready-to-use antibiotic free media for freezing and thawing of human embryos between 2PN and 4-cell stage.

- **Ready-to-use**
- **Freezing and thawing of human embryos (2PN, 4-cell stage)**
- **Contains Propandiol and Sucrose**
- **Contains Human Serum Albumin**
 Embryo Freeze - 15 g/litre
 Embryo Thaw 1 - 14 g/litre
 Embryo Thaw 2 - 14 g/litre
 Embryo Thaw 3 - 14 g/litre
- **CE marked class III (0344)**



GM501 EmbryoStore

GM501 EmbryoStore is a set of ready-to-use antibiotic free media for freezing and thawing of human embryos between 2PN and 4-cell stage.



Composition

- NaCl, KCl, KH_2PO_4 , Na_2HPO_4
- Sucrose
- Propandiol
- Human Serum Albumin

Tested specifications

- pH
- Osmolality (EmbryoStore Thaw 3)
- Sterility
- Endotoxins
- MEA

Instructions for use

Ensure all media are mixed well and warmed up to room temperature before use.

Recommended application

For further information, please see our IVF Media Recommended Usage Guide.

Order codes	Size	Storage	Shelf life*
FGY4EMF01_P_F	1 x 10 ml Freeze	2 - 8 °C	18 months
FGY4EMF01_P_T1	1 x 10 ml Thaw 1	2 - 8 °C	18 months
FGY4EMF01_P_T2	1 x 10 ml Thaw 2	2 - 8 °C	18 months
FGY4EMF01_P_T3	1 x 10 ml Thaw 3	2 - 8 °C	18 months

* from time of manufacture

Gi115/V3

Planer Limited

110 Windmill Road, Sunbury-On-Thames
Middlesex TW16 7HD, United Kingdom
A Hamilton Thorne Company

Tel: +44 (0)1932 755 000
enquiries@planer.com
www.planer.com

Manufactured by

Gynemed GmbH & Co. KG
Lübecker Straße 9 23738 Lensahn
Germany

Tel.: +49 (0) 4363 90 32 90
info@gynemed.de
www.gynemed.de