



Incubator Management Software

PIMS

Easy to use management software for Planer BT37 and CT37stax™ incubators.

Planer Incubator Management Software (PIMS) operates as both the full control interface for a single or multiple CT37stax incubators but also can be used as a management system for Planer BT37 unit.

The intuitive suite provides the convenience of controlling multiple CT37stax™ units as well monitoring key parameters of multiple BT37 units from single or multiple PCs anywhere in the laboratory.

The PIMS system is suitable for 21 CFR Part 11 compliance and fully password protected, logging actions taken by each authorised user for audit reports.

The PIMS system will alert users of any out of range alarm conditions triggered by any Planer incubator connected to the network. It allows users to acknowledge and silence alarms giving them time to visit the incubator to correct the alarm condition.

Electronic location information for each sample can be recorded via the PIMS application to assist users to easily locate specific samples and ensure they are returned back to the same location.


PIMS can graphically display real time information for each incubator connected to the system and provide full audit reports.

- Free to download
- Manage all your Planer incubators at a single point
- Real time graphical display
- Data logging of all incubator parameters
- User management
- Alarm notification management
- Suitable for 21 CFR Part 11 compliance

Want to find out more about incubation? enquiries@planer.com

Planer Incubator Management Software

PIMS




Manage all Planer incubators from a single desktop.



Display alerts and manage alarm notifications for any incubator connected to the system.



Real time graphical display of each incubator's performance.



Electronic dish labelling.

ID	Name	Admin	Active	Password
MANAGER	Manager	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	***
admin	An administrator	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	***
jdoe	John Doe	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	***
user	A User	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	***

User and password management.

Specifications may change without notice. Third party trademarks acknowledged.

BI056/V3

Planer Limited
110 Windmill Road, Sunbury-On-Thames
Middlesex TW16 7HD, United Kingdom
A Hamilton Thorne Company

Tel: +44 (0)1932 755 000
enquiries@planer.com
www.planer.com