



Cryogenic RFID Witnessing System

BLUETOOTH PROXIMITY READER

THE BLUETOOTH PROXIMITY READER IS THE NEXT ITERATION OF THE PROXIMITY READER SERIES.

The new bluetooth pepper C1 RFID reader module now provides a bluetooth interface to the PC thus significantly reducing connection dropouts and improving read times. It now also uses a USB A port to charge the internal battery.

SPECIFICATIONS

Bluetooth Interface	IEEE 802.15.1 2.4 GHz frequency Bluetooth SPP profile, BLE HID, and custom BLE service
Battery	Lithium-Ion Rechargeable Battery
Dimensions L x W x H (mm)	305x60x45
Weight	308g
LED	Green LED indicated if unit is powered on Flashing Red LED indicated if RFID tag is read
RFID	ISO15693 High Frequency 13.56 MHz I-CODE1 and ISO/IEC 15693 protocols
Data Rates	9.6 kbit/s
Protection Rating	IP 54
Operating Temperature	-40 to +80 °C
Storage Temperature	0 to +35 °C
Case Material	Nylon 12
Software	Drivers compatible with Windows 10 and 11
Warranty	1 year

WANT TO FIND OUT MORE ABOUT CRYOGENIC RFID? CONTACT US.

Specifications may change without notice. Third party trademarks acknowledged.

Si211/V1
© 2023 Planer Ltd

PLANER
PRESERVE / PROTECT / NURTURE

A Hamilton Thorne Company

www.planer.com

Manufactured by:
Cryogatt Systems Ltd
7, Abbots Business Park
Primrose Hill, Kings Langley
WD4 8FR, United Kingdom

Planer Limited
110 Windmill Road
Sunbury-On-Thames
TW16 7HD, United Kingdom
+44 (0)1932 755 000

Contact us: planer.com/contact