



Microtools

DENU-Tips

The DENU-Tips are used for the manipulation and transfer of oocytes and embryos during IVF/ ICSI procedures.

- Made out of Polyamid
- Singly and sterile packed of 20
- Sterilised for 3 years*
- Different sizes available
- Mouse-Embryo-tested
- CE marked (Class IIa)

* from time of manufacture



DENU-Tips

The DENU-Tips are used for the manipulation and transfer of oocytes and embryos during IVF/ ICSI procedures.

The DENU-Tips are ordered by supplying the Code (GY), Inner Diameter size and Number of Tips in Box as shown below.



Specification and quality control

- Considered to be Bisphenol A (BPA) free.
- The DENU-Tips are available in a variety of sizes with an inner diameter from 130 µm to 550µm. For easy differentiation, each size is color coded.
- A Mouse Embryo Test (MEA) result is available for each lot number upon request from our website.

Order codes	Inner diameter µm	Colour	Box (pieces)
FGYGY-135-20	135 µm	Yellow	20
FGYGY-140-20	140 µm	White	20
FGYGY-145-20	145 µm	White	20
FGYGY-150-20	150 µm	Green	20
FGYGY-155-20	155 µm	Green	20
FGYGY-165-20	165 µm	Black	20
FGYGY-170-20	170 µm	Red	20
FGYGY-175-20	175 µm	Red	20
FGYGY-200-20	200 µm	Orange	20
FGYGY-275-20	275 µm	Gray	20
FGYGY-300-20	300 µm	Brown	20

Gi122/V4