



## PLANER MRV



## PLANER MR7

Dimensions (H x W x D mm)	80 x 220 x 350 Mounted to power adaptor* 115 x 220 x 340	134 x 242 x 60 Mounted to desk stand 228 x 242 x 205
Weight	2.6 kg	2.4 kg
Display	240 x 64 LCD with CCFL backlight	7" TFT LCD display 800 x 480 pixels (157 x 89 mm)
Reporting	320/640 Dot thermal printer	Digital reporting and automatic filing to any network via direct Wi-Fi connection
Keypad	20 Key membrane	Resistive touchscreen
User Interface	MRV specific **	Variant of DeltaT v7 **
Programmable Cooling Rate Range	-0.01 °C/min to -99.9 °C/min	
Number of Profiles	10	Unlimited***
Steps per Profile	Up to 32	
Number of Stored Runs	5	Unlimited***
Power: Kryo 360	Powered direct from freezer	
Power: Kryo 560	Powered via separate power adaptor	DC adaptor in desk stand
Connectivity	RS232 Serial connection to PC	WiFi and Ethernet

\* When used with Kryo 560. \*\* DeltaT v7 software can be used as user interface running from a PC. \*\*\* SD Card has 8GB storage, enough for 6 million profiles or 80,000 Kryofiles.