

# THE GYNEMED MEDIALINE RANGE

## Oocyte Handling



**GM501 Flush**  
Aspiration medium



**GM501 Wash or GM501 Wash with Phenolred and Gentamicin**  
Washing medium for oocytes



**GM501 Hyaluronidase**  
Denudation of oocyte



**GM508 CultActive**  
Activation of oocytes (ionophore)

## Sperm Handling



**GM501 Gradient**  
Density gradient medium



**GM501 SpermAir or SpermActive**  
Sperm processing



**Bromelain in Dulbecco's PBS**  
Liquefy viscous semen (Bromelain)



**GM501 PVP**  
Immobilisation of sperm cells

## TESE Handling



**GM501 Collagenase**  
Digestion of testicular tissue



**GM501 SpermMobil**  
Sperm activation (Theophylline)



**SemenStain**  
Examine morphology and vitality of sperm cells



**SemenMar**  
Detect anti-sperm antibodies



**SemenLeu**  
Test elevated levels of leukocytes

## CULTURE



**GM501 Cult**  
Culture medium with Gentamicin or Gentamicin and Phenolred optional



**GM501 Mineral Oil**  
Cover of culture medium

## Cryo



**GM501 EmbryoStore**  
Slow freezing kit



**GM501 VitriStore**  
Vitrification kit



**GM501 SpermStore**  
Freezing medium (sperm and testicular tissue)