



Cryogenic RFID Witnessing System

10 X 10 BOX READER MK2

THE 10X10 BOX READER MK2 IS THE NEXT ITERATION OF THE 10X10 BOX READER SERIES.

It now uses a single USB B port which serves as both the power supply and data transfer. This makes the system more flexible as it can run off the laptop battery.

SPECIFICATIONS

Communication	USB (USB A to USB B)
Power Requirement	Power via USB
Dimensions L x W x H (mm)	180 x 180 x 40
Weight	564g
LED	Green LED indicated if unit is powered on
RFID	ISO15693 High Frequency 13.56 MHz I-CODE1 and ISO/IEC 15693 protocols
Data Rates	424 kbit/s
Protection Rating	IP 54
Operating Temperature	-40 to +80 °C
Storage Temperature	0 to +35 °C
Case Material	Nylon 12
Software	Drivers compatible with Windows 10 and 11
Warranty	1 year

WANT TO FIND OUT MORE ABOUT CRYOGENIC RFID? CONTACT US.

Specifications may change without notice. Third party trademarks acknowledged.

Si210/V1
© 2023 Planer Ltd

PLANER
PRESERVE / PROTECT / NURTURE

A Hamilton Thorne Company

www.planer.com

Manufactured by:
Cryogatt Systems Ltd
7, Abbots Business Park
Primrose Hill, Kings Langley
WD4 8FR, United Kingdom

Planer Limited
110 Windmill Road
Sunbury-On-Thames
TW16 7HD, United Kingdom
+44 (0)1932 755 000

Contact us: planer.com/contact