



Microinjectors

IM-21VM Oil Microinjector

Quick response oil-type injector

This microinjector uses oil as mediator transmitting pressure suitable for procedures such as sperm injection and holding of an oocyte.

This compact designed injector comes equipped with advanced functions that makes it easier to fill oil and remove air-bubbles.

- Simplifies troublesome oil filling and air bubble removal
- Easy tube connection
- Built-in plastic cylinder providing clear visibility



IM-21VM Oil Microinjector

This microinjector uses oil as mediator transmitting pressure suitable for procedures such as sperm injection and holding of an oocyte.



Order code	
FNHIM-21	IM-21VM oil microinjector

Specification	
Accessories included	IM-H3 Injection holder set Cleaning rod HI01PK01 Silicone rubber gasket IMO4 O ring Silicone grease
Movement range	Plunger 20 mm, Full rotation of knob 500 μ m
Control capacity	Full rotation of knob approx. 10 μ l
Dimensions/Weight	127 – 147 (W) \times 56 (D) \times 78 (H) mm, 620 g
Tube length	1 m

Ni191/V1

Planer Limited

110 Windmill Road, Sunbury-On-Thames
Middlesex TW16 7HD, United Kingdom
A Hamilton Thorne Company

Tel: +44 (0)1932 755 000
enquiries@planer.com
www.planer.com

Manufactured by

Narishige Scientific Instrument Lab.
27-9, Minamikarasuyama 4-chome,
Setagaya-ku, Tokyo 157-0062
Japan

Tel: +81-(0)3-3308-8233
sales@narishige.co.jp
www.narishige-group.com