



IVF Needles

Intrauterin Insemination Catheters (IUI)

IUI catheter range is designed to improve your success rate in artificial insemination.

To allow for different users' preferences, the Kitazato IUI catheter range is available in two different lengths (10 & 18 cm), two different gauges (5.2 and 6 Fr) and four different levels of rigidity (hard, ultimate, intermediate and soft).

The smooth and rounded tip allows for easy insertion and the double lateral opening tip improves the dispersion of the semen.

- Improved semen dispersion
- Easy insertion
- Large range of sizes
- Variety of catheter rigidity options



Intrauterin Insemination Catheters (IUI)

Improve your artificial insemination success rates with the Kitazato range of IUI catheters, which are available in 4 levels of rigidity and two different IUI catheter lengths.



Order Number	Description	Diameter	Length	Total volume minus hub
FDM131181	Hard - Long Type	6 Fr	18 cm	0.03 ml
FDM131186	Ultimate - Long Type	6 Fr	18 cm	0.03 ml
FDM131182	Intermediate - Long Type	6 Fr	18 cm	0.1 ml
FDM131184	Soft - Long Type	5.2 Fr	18 cm	0.1 ml
FDM140198	Stylet Cannula for Soft and Intermediate		18 cm	
FDM131101	Hard - Short Type	6 Fr	10 cm	0.015 ml
FDM131106	Ultimate - Short Type	6 Fr	10 cm	0.015 ml
FDM131102	Intermediate - Short Type	6 Fr	10 cm	0.05 ml
FDM131104	Soft - Short Type	5.2 Fr	10 cm	0.05 ml
FDM140118	Stylet Cannula for Soft and Intermediate		10 cm	

*Kitazato IUI Catheters are sold in boxes of 10 units.

Quality Control

- MEA (Mouse Embryo Assay): One cell assay \geq 80% after 96 hours
- Endotoxin: < 20EU/device
- SAL 10^{-6}
- Sterility Test (Bacteria, Fungi)
- Cytotoxicity Test
- Intracutaneous Reactivity Test
- Sensitisation Text

Composition

- Polyvinylchloride DEHP free
- Stainless Steel SUS304

Ki145/V4

PLANER
PRESERVE / PROTECT / NURTURE

A Hamilton Thorne Company

www.planer.com

Manufactured by:
Kitazato Corporation
100-10 Yanagishima
Fuji, Shizuoka
416-0932 Japan
www.kitazato-ivf.com

Planer Limited
110 Windmill Road
Sunbury-On-Thames
TW16 7HD, United Kingdom
+44 (0)1932 755 000

Contact us: planer.com/contact