



Microinjectors

IM-HD1TVM Pneumatic Injector (for Holding Side)

Injector for the holding side exclusively for TAKANOME™

This Narishige pneumatic injector is intended exclusively for holding an egg. It can be attached to the joystick of TAKANOME™, allowing the holding and releasing an egg in a delicate and secure manner operated only with a forefinger.

The pneumatic injector enables the micromanipulation and suctioning work at the same time. It can be used with one of two types of tubing, with different inner diameters, which allows the response speed to be customised.

- **Designed specifically for holding oocytes**
- **Precise hold and release, using only forefinger**
- **Response speed can be adjusted by changing tubing**



IM-HD1TVM Pneumatic Injector (for Holding Side)

This Narishige Pneumatic Injector is designed specifically for holding the oocyte.



Order code	
FNHIM-HD1T	IM-HD1 TVM pneumatic injector (for holding side of TAKONOME)

Specification	
Accessories included	Injection Holder (FNHHI-9) Cleaning Rod Tubing (FNHCT-4) Tubing (FNHCT-1-1.5) Silicone Rubber Gasket (FNHHI01PK01) Hex Wrench
Dimensions/Weight	167-214 mm (W) x 55 mm (D) x 78 mm (H) / 680 g

Ni190/V1

Planer Limited

110 Windmill Road, Sunbury-On-Thames
Middlesex TW16 7HD, United Kingdom
A Hamilton Thorne Company

Tel: +44 (0)1932 755 000
enquiries@planer.com
www.planer.com

Manufactured by

Narishige Scientific Instrument Lab.
27-9, Minamikarasuyama 4-chome,
Setagaya-ku, Tokyo 157-0062
Japan

Tel: +81-(0)3-3308-8233
sales@narishige.co.jp
www.narishige-group.com