



## Sample tracking software

### Significantly reduce the amount of time to recover stored samples with the Pro-Curo inventory management solution.

Pro-Curo is simple to get up and running so no intense training will be required. And if you do have a problem, there is around the clock telephone support.

The software can be quickly adapted to your needs – whether you need to track your specimens across multiple sites or need to add a patient consent form to sample information.

Want to find out more about tracking your samples?  
[enquiries@planer.com](mailto:enquiries@planer.com)

The software is fully FDA CFR 21 part 11 Compliant, giving you a complete history of every transaction and an auditable trail of your samples.

Sign up for a 14 day free trial to see just how easy Pro-Curo software is to get up and running to track samples in your laboratory.

- Fully scalable - with support for one lab up to multiple sites
- Convenient standard integration with bar code readers
- Build context around sample information - option to include attachments
- Stay informed with granular traceability for all samples
- Secure multi-level defined access control



# Sample tracking software

Pro-Curo



FEATURES	
Max no of samples	100,000,000
Complete history of every transaction	Yes
Multi PC	Yes - 1, 3, 5, Unlimited
Database structure	MSSQL & MySQL
User manager	Yes
Sample viewer restricted by project	Yes
Auto logoff function	Yes
Split samples function	Yes
Customisable user fields	Yes
Edit multiple samples	Yes
Book in and book out function	Yes
Multi level location structure	Yes - 7
Grid view	Yes
Standard and custom report generation	Yes
Label design	Graphical
Data import function	Yes
Sample search function	Yes
Advanced SQL search builder	Yes
Date alerts	Yes
Double blind auditing function	Yes

Specifications may change without notice, third party trademarks acknowledged. Di044/V2

© 2019 Planer Ltd

## Planer Limited

110 Windmill Road, Sunbury-On-Thames  
Middlesex TW16 7HD, United Kingdom  
A Hamilton Thorne Company

**Tel:** +44 (0)1932 755 000  
**Fax:** +44 (0)1932 755 001

enquiries@planer.com  
www.planer.com