



GYNEMED

www.gynemed.de

Content

- 1 Vial with 2 ml GM501 SpermAir
- 2 Syringes
- 3 Injection cannulae
- 1 IUI-Catheter
- Instructions for use

Insemination Kit

Insemination

Insemination Kit

The Insemination Kit is a complete system for the simple and safe preparation and processing of human spermatozoa out of the ejaculate, for homologous and heterologous intrauterine inseminations (IUI).

The Insemination Kit uses the self motility of the male germ cells to isolate motile spermatozoa in high concentrations. We recommend the use if the ejaculate is normozoospermic or slightly oligoand/ or asthenozoospermic.

- **Ready-to-use**
- **Complete system for IUI**
- **Contains GM501 SpermAir**



Insemination Kit

The Insemination Kit is a complete system for the simple and safe preparation and processing of human spermatozoa out of the ejaculate, for homologous and heterologous intrauterine inseminations (IUI).

The Insemination Kit uses the self motility of the male germ cells to isolate motile spermatozoa in high concentrations. We recommend the use if the ejaculate is normozoospermic or slightly oligoand/ or asthenozoospermic.




One kit contains

- 1 x Vial with 2 ml of GM501 SpermAir medium (sterile)
- 2 x 2 ml syringes (sterile)
- 1 x short cannula (sterile)
- 2 x long cannulas (sterile)
- 1 x Insemination catheter standard or memo (sterile)
- 1 x Instructions for use
- 1 x Ampulla rack

Recommended application

For further information, please see our IVF Media Recommended Usage Guide.

Instructions for use



Order code	Size	Storage	Shelf life*
FGY4SA-Kit-002	1 x Standard kit	2 - 8 °C	6 months

* from time of manufacture

Gi118/V3

Planer Limited

110 Windmill Road, Sunbury-On-Thames
Middlesex TW16 7HD, United Kingdom
A Hamilton Thorne Company

Tel: +44 (0)1932 755 000
enquiries@planer.com
www.planer.com

Manufactured by

Gynemed GmbH & Co. KG
Lübecker Straße 9 23738 Lensahn
Germany

Tel.: +49 (0) 4363 90 32 90
info@gynemed.de
www.gynemed.de